



# Deuterated Gangliosides for Mass Spectrometry



Deuterated gangliosides are ideal for the identification and quantification of gangliosides in samples and biological systems using mass spectrometry. In addition to Matreya's collection of D<sub>3</sub>-labeled gangliosides, our scientists have synthesized a new set of high-purity D<sub>9</sub>-labeled gangliosides that enable improved separation of the labeled *m/z* from the naturally occurring material and its abundant isotopes. With a better means for separation, this improves the quantitation of these compounds in biological samples by LC-MS/MS methods.

## Advantages of D<sub>9</sub>-Labeled Gangliosides

Natural gangliosides have a large molecular weight and a high <sup>13</sup>C isotopic abundance up to the M+5 isotope. For example, the GM<sub>1</sub> M+3 isotope is very significant but has approximately the same mass as its D<sub>3</sub>-labeled standard. This natural M+3 isotope can contribute to variability in the deuterated internal standard response as GM<sub>1</sub> levels vary from sample to sample and can negatively impact quantification of GM<sub>1</sub>. A D<sub>9</sub>-labeled standard overcomes this problem because the natural abundance of GM<sub>1</sub> with nine <sup>13</sup>C atoms is effectively zero. That is, the natural isotopes will not interfere with the GM<sub>1</sub>-D<sub>9</sub> signal.

View a comparison of the mass spectra of C16-D<sub>0</sub>-GM<sub>1</sub>, C16-D<sub>3</sub>-GM<sub>1</sub>, and C16-D<sub>9</sub>-GM<sub>1</sub> on the next page

## D<sub>9</sub>-Labeled Gangliosides

Catalog No.	Product Name	Formula Weight	Purity
2057	N-Hexadecanoyl-D <sub>9</sub> (13,13,14,14,15,15,16,16,16)-monosialoganglioside GM <sub>1</sub> (NH <sub>4</sub> <sup>+</sup> salt)	1528 + NH <sub>3</sub>	98+%
2058	N-Hexadecanoyl-D <sub>9</sub> (13,13,14,14,15,15,16,16,16)-monosialoganglioside GM <sub>2</sub> (NH <sub>4</sub> <sup>+</sup> salt)	1366 + NH <sub>3</sub>	98+%
2059	N-Hexadecanoyl-D <sub>9</sub> (13,13,14,14,15,15,16,16,16)-monosialoganglioside GM <sub>3</sub> (NH <sub>4</sub> <sup>+</sup> salt)	1163 + NH <sub>3</sub>	98+%

## D<sub>3</sub>-Labeled Gangliosides

Catalog No.	Product Name	Formula Weight	Purity
2050	N-omega-CD <sub>3</sub> -Octadecanoyl monosialoganglioside GM <sub>1</sub> (NH <sub>4</sub> <sup>+</sup> salt)	1550 + NH <sub>3</sub>	98+%
2051	N-omega-CD <sub>3</sub> -Octadecanoyl monosialoganglioside GM <sub>2</sub> (NH <sub>4</sub> <sup>+</sup> salt)	1388 + NH <sub>3</sub>	98+%
2052	N-omega-CD <sub>3</sub> -Octadecanoyl monosialoganglioside GM <sub>3</sub> (NH <sub>4</sub> <sup>+</sup> salt)	1185 + NH <sub>3</sub>	98+%
2054	N-omega-CD <sub>3</sub> -Octadecanoyl disialoganglioside GD <sub>3</sub>	1476	98+%

MATREYA, LLC  
WWW.MATREYA.COM


Available in Canada from...

MJS  
**BioLynx**  
INC.


1-888-593-5969 • www.biolyx.ca • tech@biolyx.ca




### N-Hexadecanoyl monosialoganglioside GM<sub>1</sub>



### N-*omega*-CD<sub>3</sub>-Hexadecanoyl monosialoganglioside GM<sub>1</sub>



### N-Hexadecanoyl-D<sub>9</sub> (13,13,14,14,15,15,16,16,16)- monosialoganglioside GM<sub>1</sub>



Natural isotopic abundance of GM<sub>1</sub> (top spectrum) compared to the D<sub>3</sub>-labeled (middle spectrum) and D<sub>9</sub>-labeled (bottom spectrum) versions.

 MATREYA, LLC  
WWW.MATREYA.COM

If you cannot find a particular Deuterated Ganglioside,  
please contact our Technical Team for a custom synthesis estimate.

MJS  
BioLynx  
INC.

1-888-593-5969 • [www.biolynx.ca](http://www.biolynx.ca) • [tech@biolynx.ca](mailto:tech@biolynx.ca)