

Development of a high-throughput CUT&RUN platform for epigenomic mapping of rare primary immune cells

Keith Maier¹, Ellen N. Weinzapfel¹, Allison R. Hickman¹, Matthew R. Marunde¹, Carolina P. Lin¹, Danielle N. Maryanski¹, Liz Albertorio-Saez¹, Dughan J. Ahimovic², Michael Bale², Juliana J. Lee³, Michael-Christopher Keogh¹, Immunological Genome Project Consortium

¹EpiCypher Inc., Research Triangle Park, Durham NC 27709, USA. ²Laboratory of Epigenetics and Immunity, Department of Pathology and Laboratory Medicine, Weill Cornell Medicine, New York, NY, USA. ³Department of Immunology, Harvard Medical School, Boston, MA, USA.



Epigenetic regulation is central to cell and gene therapy, but has been challenging to study

- Many genomic strategies for cell & gene therapy focus on transcription; however, RNA-seq reveals the **outcomes** – not driving **mechanisms**
- Epigenomics is the solution:** Mapping the location of histone post-translational modifications (PTMs) and chromatin-associated proteins, such as transcription factors, provides molecular insights that are central to cell fate and function
- However, existing epigenomic technologies, such as **ChIP-seq**, are limited by high costs, poor sensitivity & reliability, and complicated sample prep
- These challenges have precluded epigenomic analysis for cell & gene therapy




Figure 1: Understanding epigenetic regulation is critical to successful cell and gene therapy applications:

- iPSCs
- CAR T-cells
- T cell exhaustion
- dCas9/Cas9 targeting

CUTANA™ CUT&RUN assays provide important advantages compared to ChIP-seq





Figure 2: Overview of the CUTANA CUT&RUN workflow and advantages compared to ChIP-seq. Because CUT&RUN releases antibody-bound fragments into solution (A), it has improved signal-to-noise even with significantly reduced cell numbers and sequencing depth (B).

Automated CUTANA CUT&RUN Assays: high-throughput, customized epigenomics



autoCUT&RUN enables ultra-sensitive epigenomic profiling from low cell numbers





Figure 4: Automated CUTANA CUT&RUN (autoCUT&RUN) was used to generate maps for various histone PTMs using decreasing amounts of K562 cells (A). Pearson correlation matrix of H3K4me3 data shows high concordance across cell numbers for each target (B).

The key to autoCUT&RUN: SNAP-CUTANA™ Spike-in Controls and highly specific antibodies



Unique validation strategy enables reliable profiling of transcription factors and chromatin-associated proteins





Figure 6: Our stringent validation protocol involves multiple genome-wide analyses to ensure: Peak and motif enrichment consistent with biological function (A, C), High signal over background in genome-wide analyses (B), Reliable performance across production lots, FOXA1 antibody (EpiCypher 13-2001) data are shown

High quality antibodies enable low-input autoCUT&RUN experiments



Proof-of-concept: Robust profiling using 10,000 mouse immune cells




Figure 8: autoCUT&RUN is reliable for immune cell profiling. autoCUT&RUN was performed using 10,000 lightly cross-linked mouse B cell nuclei (provided by the ImmGen Consortium) and same antibodies as in Figure 7.

Conclusions

- Epigenomics is central to understanding gene regulatory processes, but historical methods (ChIP-seq) are unreliable
- Ultra-sensitive CUTANA CUT&RUN assays are poised to dramatically change the field, improving access to high-resolution chromatin mapping
- EpiCypher developed **automated CUT&RUN for high throughput and cost-effective chromatin mapping**, ideal for studying cell & gene therapy at scale
- These efforts are bolstered by **quantitative SNAP-CUTANA Spike-in Controls and highly specific antibodies**, both of which were crucial to optimizing autoCUT&RUN for immune cell profiling

autoCUT&RUN defines immune cell differentiation pathways for advanced cell & gene therapy research

- Collaboration with ImmGen Consortium
- Each assay only required 10,000 FACS-sorted cells
- Compatible with stimulated and cross-linked cells
- Standardized sample handling improved yields
- Cell-type specific peaks detected – useful for studying cell differentiation

Uncover cell type-specific gene regulatory profiles




Figure 9: autoCUT&RUN reveals distinct H3K4me3 (marks active promoters) and H3K27me3 (denotes repressed genes) profiles across FACS-sorted primary mouse granulocytes, type 3 innate lymphoid cells, and nature killer cells (Ly49H+), provided by ImmGen. 10,000 nuclei were used per reaction.

In-depth studies of individual immune cell populations reveals novel biology




Figure 10: High-resolution profiling of FACS-sorted type 3 ILCs using autoCUT&RUN identifies unique genomic compartments, including active regulatory elements (H3K4me1, H3K27ac), promoters (H3K4me3), and gene bodies (H3K36me3), as well as repressed genes (H3K27me3) and transcription factor binding sites (CTCF). 10,000 nuclei (from ImmGen) were used per reaction.

Applications of CUTANA assays

- Exhausted T cells (PMID: 35930654)
- CAR T cell expansion (PMID: 36944333)
- dCas9/Cas9 targeting (PMID: 35849129)



CHROMATOGRAPHIC SPECIALTIES INC.

1-888-593-5969 • biolynx.ca • tech@biolynx.ca